**$\alpha$ -AMYLASE-SL ASSAY**

**CATALOGUE NUMBER:** 341-10      **SIZE:** 10 x 10 mL  
 341-40      4 x 100 mL

**INTENDED USE**

For the IN VITRO quantitative measurement of  $\alpha$ -amylase in serum.

**TEST SUMMARY**

$\alpha$ -Amylase (1,4- $\alpha$ -D-Glucan glucanohydrolase, EC 3.2.1.1) is an enzyme of the digestive tract. It hydrolyzes dietary starch and glycogen to form maltose by splitting their chains at alternate hemiacetal linkages. The enzyme is normally secreted into the digestive tract from the parotid glands and the pancreas. In diseases affecting these glands, and particularly when the pancreatic duct is obstructed, the amount of enzyme in serum is increased.<sup>(1)</sup>

Methods employed for assaying amylase activity are:

1. Saccharogenic assays, where the amount of reducing sugars formed are measured after serum and a starch substrate are incubated.
2. Amyloclastic methods, which follow the decrease in concentration of the starch substrate when it is incubated with serum.
3. Enzymatic methods in which the maltose formed is metabolized (usually through glucose) in a series of reactions which generate NADH.
- 4a. Enzymatic methods which quantitate the release of p-nitrophenol from the p-nitrophenol maltopentaoside and hexaoside substrates.<sup>(2)</sup>
- 4b. Enzymatic method which quantitates the release of 2-chloro-p-nitrophenol from the 2-chloro-p-nitrophenyl- $\alpha$ -D-maltotrioside (CNP3) substrate.<sup>(3)</sup>

This amylase assay is a one step enzymatic procedure which uses the substrate CNPG3.<sup>(3)</sup> The advantage being that the  $\alpha$ -amylase reacts directly, without auxiliary enzymes, with the CNPG3 substrate releasing 2-chloro-p-nitrophenol.

**TEST PRINCIPLE**

$\alpha$ -Amylase

CNP3 (2-chloro-p-nitrophenyl- $\alpha$ -D-maltotrioside)  $\rightarrow$  2-chloro-p-nitrophenol

The  $\alpha$ -amylase reacts directly with CNPG3 to release 2-chloro-p-nitrophenol which is monitored at 405 nm. The rate of increase in absorbance at 405 nm is proportional to the  $\alpha$ -amylase in the sample.

**REAGENTS**

$\alpha$ -Amylase Reagent: A solution containing at least 2.25 mmol/L CNPG3, buffer (pH 6.0 at 25°C), activators, and stabilizers.

**WARNINGS & PRECAUTIONS FOR USE**

S24/25: Avoid ingestion and contact with skin and eyes.  
 See Material Safety Data Sheet for additional information.  
 Saliva contains amylase. Avoid contamination of testing materials or containers with saliva.

**REAGENT PREPARATION, STORAGE, AND STABILITY**

The reagent is provided in a ready to use format.

The reagent included is stable until the expiry date stated on the vial label when stored at 2-8°C.

Stability claims are based on real time studies.

**REAGENT DETERIORATION**

The reagent solution should be clear. Turbidity would indicate deterioration.

**DISPOSAL**

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

**SPECIMEN**

Fresh, clear, unhemolysed serum.

**SAMPLE STORAGE**

Amylase is stable in serum for at least one week when stored at room temperature or for several months at 4°C.<sup>(4)</sup>

**ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(6)</sup>**

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Amylase in saliva will falsely elevate test results if testing materials become contaminated with saliva.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
121 U/L	Hemoglobin	200 mg/dL	31 $\mu$ mol/L
121 U/L	Bilirubin	40 mg/dL	684 $\mu$ mol/L
121 U/L	Intralipid	600 mg/dL	1800 mg/dL (20.3 mmol/L) Simulated Triglycerides

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(5)</sup>

**ANALYTICAL PROCEDURE****MATERIALS PROVIDED**

Genzyme Diagnostics'  $\alpha$ -Amylase reagent.

**MATERIALS REQUIRED (BUT NOT PROVIDED)**

- 1) Automated analyzer capable of accurately measuring absorbances at appropriate wavelength as per instrument application.
- 2) Calibration material (if applicable).
- 3) Quality Control materials.

**TEST CONDITION**

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:60, and wavelength reading of (primary/ secondary) 416 nm / 660 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

**CALIBRATION**

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

**QUALITY CONTROL**

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

**CALCULATIONS**

The analyzer automatically calculates the  $\alpha$ -amylase concentration of each sample.

**TEST LIMITATIONS**

A sample with an  $\alpha$ -amylase concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

**REFERENCE INTERVALS<sup>(4)</sup>**

The expected values of  $\alpha$ -amylase were determined using serum obtained from 120 apparently normal adults following the procedures outlined in CLSI Document C28-P. The range was determined to be 18-87 U/L at 37°C.

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

**PERFORMANCE CHARACTERISTICS**

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

**RESULTS**

$\alpha$ -Amylase activity is reported as U/L.

## REPORTABLE RANGE (CLSI EP6)<sup>(6)</sup>

The linearity of the procedure described is 2000 U/L. The lower limit of detection of the procedure is 3 U/L. The data results in a reportable range of 3-2000 U/L.

## PRECISION STUDIES (CLSI EP5)<sup>(6)</sup>

Data was collected on two concentrations of control sera using a single lot of reagent in 40 runs conducted over 20 days.

Concentration (U/L)	Total SD (U/L)	Total CV (%)	Within-Run SD (U/L)	Within-Run CV (%)
46	1.5	3.3	1.2	2.7
390	8.8	2.3	2.5	0.6

## ACCURACY (CLSI EP9)<sup>(6)</sup>

The performance of this method (y) was compared with the performance of another CNPG3 method (x) on a Roche/Hitachi® analyzer. Forty samples ranging from 16-1779 U/L gave a correlation coefficient of 0.9984. Linear regression analysis gave the following equation:

$$\text{This method} = 0.9603 (\text{reference method}) + 4.5 \text{ U/L.}$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

### The Americas

Genzyme Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: [customerservice@genzyme.com](mailto:customerservice@genzyme.com)

[peidiagnosticttechnical@genzyme.com](mailto:peidiagnosticttechnical@genzyme.com)

### International

Genzyme Diagnostics  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

Email: [ukdiagcustomerservice@genzyme.com](mailto:ukdiagcustomerservice@genzyme.com)

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