## BILIRUBIN-D ASSAY

**CATALOGUE NUMBER:** 202-S7 **SIZE:** 2 x 125 mL, 1 x 7.5 mL

### INTENDED USE

For the IN VITRO quantitative measurement of bilirubin (direct) concentration in serum.

### TEST SUMMARY

Bilirubin is a bile pigment normally found in serum as a result of red cell destruction. It is a product of hemoglobin breakdown by the reticuloendothelial system and exists in two forms. Unconjugated (indirect) bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted.

Elevated serum bilirubin is indicative of impairment of the liver, excessive hemolysis, or obstruction of the biliary tract.

Ehrlich<sup>(1)</sup> described the coupling of bilirubin with diazotized sulfanilic acid to form a red pigment in neutral solutions and a blue pigment in strongly acid or alkaline solutions. Conjugated (direct-reacting) bilirubin couples directly with the reagent in an aqueous medium. An accelerator can be added to promote the coupling of free bilirubin with the diazo reagent. This allows the measurement of total bilirubin. Alcohol<sup>(2)</sup> or caffeine-sodium benzoate<sup>(3)</sup> are commonly used as accelerators. Winsten and Cehelyk<sup>(4)</sup> used dimethylsulfoxide (DMSO) as an accelerator. This method is based on the Walters and Gerarde<sup>(5)</sup> modification of the DMSO method.

### TEST PRINCIPLE

Direct Bilirubin + Diazotized Sulfanilic Acid → Azobilirubin (Red)

Sulfanilic acid reacts with sodium nitrite to form diazotized sulfanilic acid. Direct bilirubin reacts with diazotized sulfanilic acid forming azobilirubin.

The absorbance of the reaction mixture at 555 nm is directly proportional to the concentration of direct bilirubin.

### REAGENTS

Bilirubin Reagent (Direct) (R1): A solution containing 35.6 mmol/L sulfanilic acid and hydrochloric acid.

Sodium Nitrite Reagent (R2): A solution containing 43.5 mmol/L sodium nitrite.

### WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.

Avoid ingestion.

The reagent is corrosive. Use gloves when handling.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

A direct bilirubin working reagent may be prepared by adding one drop (0.05 mL) of nitrite reagent to each 3.0 mL of direct bilirubin reagent.

Supplied reagents are stable at 18-26°C until expiry date.

Working reagent is stable for 8 hours when stored in an amber bottle at 18-26°C.

Reagents should be kept out of direct sunlight. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration. The nitrite reagent should be discarded if a dark yellow discoloration is present.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### SAMPLE STORAGE

The specimen must be protected from light. Samples should be stored at 2-8°C and are stable for seven days. Samples may be stored at minus 20-0°C for three months.<sup>(6)</sup>

### ANALYTICAL SPECIFICITY

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

A hemoglobin concentration of greater than 5 g/L in the sample can interfere with the analysis.

A summary of the influence of drugs on clinical laboratory test may be found by consulting Young, D.S.<sup>(7)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Genzyme Diagnostics' Bilirubin (direct) reagents.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Spectrophotometer or analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material
- 3) Quality Control materials

#### TEST CONDITION

For data presented in this insert, studies using this reagent were performed using an endpoint test mode, with a sample to reagent ratio of 1:15 and a wavelength reading of 555 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the bilirubin (direct) concentration of each sample.

#### TEST LIMITATIONS

A sample with a bilirubin (direct) concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(6)</sup>

0.0-0.2 mg/dL (0.0-4.0 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

These performance characteristics were generated in Genzyme Diagnostics laboratories.

#### RESULTS

Bilirubin (direct) concentration is reported as mg/dL (µmol/L).

#### REPORTABLE RANGE

The linearity of the procedure described is 20 mg/dL (342 µmol/L). Linearity using automated procedures will depend on the sample/reagent ratio used.

## PRECISION STUDIES

Total precision data was collected on two concentrations of control sera in 20 runs conducted over 10 days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	µmol/L	mg/dL	µmol/L		mg/dL	µmol/L	mg/dL	µmol/L	
2.5	43.4	0.05	0.9	2.1	0.4	7.0	0.03	0.5	7.1
15.5	264.5	0.69	11.8	4.5	2.5	42.8	0.04	0.7	1.6

Within run precision data was collected on two concentrations of control sera each run 54 times in a single assay.

## ACCURACY

The performance of this method (y) was compared with the performance of another commercially available method (x) using 54 patient serum samples. The correlation coefficient was 0.99. Linear regression analysis gave the following equation:

$$\text{This method} = 1.02 (\text{reference method}) - 0.0 \text{ mg/dL} (0.04 \text{ µmol/L})$$

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

## TRADEMARKS

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**The Americas**  
Genzyme Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: customerservice@genzyme.com  
peidiagnosticttechnical@genzyme.com

**International**  
Genzyme Diagnostics  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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