ETHANOL ASSAY

CATALOGUE NUMBER: 229-29 **SIZE:** 12 x 2.3 mL + 3 x 125 mL + 1 x 15 mL

INTENDED USE


For the IN VITRO quantitative measurement of ethanol in serum and plasma.

TEST SUMMARY

Ethanol (ethyl alcohol) is the most common toxic substance involved in medical-legal cases. Ethanol consumption is often a factor in all types of accidents and ethanol poisoning can be fatal. Ethanol testing of comatose patients aids in the differential diagnosis of the cause of the coma.⁽¹⁾

Early methods for ethanol analysis depended on separation of the ethanol from the specimen. The method described here is an enzymatic procedure which does not require separation and is rapid and easily automated. This method is based on the first enzymatic method described by Bonnichsen and Theorelle⁽²⁾ and later modified by other authors.⁽³⁾

TEST PRINCIPLE

The enzyme alcohol dehydrogenase (ADH) catalyzes the conversion of ethanol to acetaldehyde and NADH. The formation of acetaldehyde is favoured when the reaction takes place at pH 9. To ensure that the reaction goes to completion, the acetaldehyde formed is removed from the system by its reaction with hydrazine sulfate. The increase in absorbance at 340 nm due to the reduction of NAD to NADH is proportional to the amount of ethanol present.

REAGENTS

NAD-ADH Reagent: A solution (after reconstitution) containing 8.69 mmol/L NAD, 1,450 U/mL alcohol dehydrogenase (yeast), and 9.9 mmol/L glutathione.

Ethanol Buffer: A solution containing 75 mmol/L tetrasodium pyrophosphate, 22 mmol/L glycine, and 75 mmol/L hydrazine sulfate.

Ethanol Calibrator: A solution containing 92 mg/dL or 0.092% w/v (20 mmol/L) ethanol.

WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.
Avoid ingestion.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE, AND STABILITY

Add 2.3 mL of deionized water to the required number of vials of NAD-ADH reagent. Allow one minute for reconstitution, then mix gently.

A single working reagent may be prepared by mixing the reconstituted NAD-ADH reagent and the ethanol buffer in a ratio of 1:10.

Supplied reagents are stable at 2-8°C until expiry date.

The reconstituted NAD-ADH reagent is stable at 2-8°C for seven days or for two months at minus 20-0°C. The working reagent is stable at 2-8°C for seven days.

The ethanol calibrator should be tightly sealed when not in use, as exposure to the atmosphere causes evaporation of the ethanol.

Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum or plasma. Specimens should be kept sealed.

SAMPLE STORAGE

Samples may be stored refrigerated for up to 2 weeks.⁽⁴⁾ Sample should be at room temperature (18-26°C) for testing. Always tightly cap sample tubes to prevent evaporation of alcohol.

Ethanol should not be used to swab the venipuncture site or to clean or sterilize glassware or other equipment used to collect the specimen or perform the assay. Avoid the use of test tubes, cuvettes, or pipettes which have been cleaned or sterilized with ethanol.

ANALYTICAL SPECIFICITY

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Testing was performed in Genzyme Diagnostics laboratories to determine the specificity of the ethanol reagent. Solutions of potentially interfering compounds were prepared and assayed. The concentrations of the solvents and the assay results, expressed as units of ethanol, are listed below.

| Solvent | Concentration (w/v) | Units of Ethanol | | | % Reactivity |
|-----------------|---------------------|------------------|-------|--------|--------------|
| | | % w/v | mg/dL | mmol/L | |
| Acetone | 0.1 | 0.000 | 0 | 0.0 | 0% |
| | 2.0 | 0.000 | 0 | 0.0 | 0% |
| n-Butanol | 0.1 | 0.025 | 25 | 5.4 | 25% |
| | 2.0 | 0.311 | 311 | 67.5 | 16% |
| Ethylene Glycol | 0.1 | 0.001 | 1 | 0.2 | 1% |
| | 2.0 | 0.014 | 14 | 3.0 | 0.7% |
| Methanol | 0.1 | 0.000 | 0 | 0.0 | 0% |
| | 2.0 | 0.001 | 1 | 0.2 | 0.05% |
| Isopropanol | 0.1 | 0.005 | 5 | 1.1 | 5% |
| | 2.0 | 0.085 | 85 | 18.5 | 4% |

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Genzyme Diagnostics' Ethanol reagents and calibrator.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Quality Control materials.
- 3) Deionized water for reconstitution

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:123, and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S, please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

An ethanol calibrator is included and should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the ethanol concentration of each sample.

TEST LIMITATIONS

A sample with an ethanol concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS

Values up to 10 mg/dL (0.01% w/v or 2.2 mmol/L) may be observed in sera due to the presence of volatile reducing substances and do not necessarily indicate ingestion of ethanol.⁽⁶⁾

The legal limit for drivers is less than 80 mg/dL (0.08% w/v or 17 mmol/L).⁽⁷⁾

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

RESULTS

Ethanol concentration is reported as mg/dL (% w/v or mmol/L).

REPORTABLE RANGE

The linearity of the procedure described is 230 mg/dL (0.23% or 50 mmol/L).

PRECISION STUDIES

Day to day precision was established by assaying three spiked pooled sera twice a day for 10 days.

| Concentration | | Total SD | | Total CV % | Concentration | | Within Run SD | | Within Run CV % |
|---------------|--------|----------|--------|------------|---------------|--------|---------------|--------|-----------------|
| mg/dL | mmol/L | mg/dL | mmol/L | | mg/dL | mmol/L | mg/dL | mmol/L | |
| 52 | 11.3 | 0.8 | 0.17 | 1.5 | 50 | 10.9 | 0.3 | 0.07 | 0.6 |
| 97 | 21.1 | 2.2 | 0.48 | 2.3 | 99 | 21.5 | 1.1 | 0.24 | 1.1 |
| 137 | 29.7 | 3.4 | 0.74 | 2.5 | 150 | 32.6 | 7.0 | 1.52 | 4.7 |

Within day precision was established by assaying three spiked sera twenty times.

ACCURACY (Serum)

The performance of this method (y) was compared with the performance of a similar UV alcohol dehydrogenase method (x). Linear regression analysis gave the following equation:

$$\text{This method} = 0.991 (\text{reference method}) + 1.2 \text{ mg/dL or } 0.0012\% \text{ w/v (0.26 mmol/L)}.$$

ACCURACY (Plasma)

The performance of this method (y) was compared with the performance of a similar UV alcohol dehydrogenase method (x). Linear regression analysis gave the following equation:

$$\text{This method} = 0.9984 (\text{reference method}) - 1.2 \text{ mg/dL or } 0.0012\% \text{ w/v (0.26 mmol/L)}.$$

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

TRADEMARKS

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

Diagnostics

The Americas

Genzyme Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: customerservice@genzyme.com

peidiagnosticttechnical@genzyme.com

www.genzymediagnosics.com

International

Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Internacional
Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, RU

Teléfono: 800-565-0265
Fax: 902-628-6504
Correo electrónico:
customerservice@genzyme.com
pediagnosticttechnical@genzyme.com

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Burtis, C.A., Ashwood, E.R. (Ed.), Tietz, Textbook of Clinical Chemistry, W.B. Saunders Co., Toronto, 1171-1175 (1994).
2. Bonnichsen, R.K., Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
3. Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues. Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
4. Tietz, N.W., (Editor), Clinical Guide to Laboratory Tests, 3rd Edition, 1995.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, 1990.
6. Tietz, N.W., (Editor), Fundamentals of Clinical Chemistry, W.B. Saunders Company, Philadelphia, 1970, p. 844.
7. McQueen, M.J., (Editor), SI Units in Clinical Chemistry (1980), Clinical Chemistry section Health Sciences Laboratory Medicine Program.

IN22929-8
January 28, 2010