

OSOM[®]

Influenza A&B Test



FOR USE WITH NASAL SWAB SPECIMENS


FEATURES:

- ✓ **ACCURATE, DIFFERENTIAL DIAGNOSIS;** Differentiates between influenza A&B
- ✓ **QC INSIDE;** Includes 2 additional test sticks for external QC testing
- ✓ **NASAL SWAB SPECIMENS**
- ✓ **FAST;** Less than 1 minute hands on time; results in 10 minutes or less
- ✓ **ROOM TEMPERATURE STORAGE**
- ✓ **DETECTS H1N1 STRAINS***

*Although this test has been shown to detect the Mexico/4108/2009 H1N1 virus cultured from a positive human respiratory specimen, the performance characteristics of this device with clinical specimens that are positive for the 2009 H1N1 influenza virus have not been established.



TEST PROCEDURE / RESULTS



ORDERING INFORMATION

- OSOM[®] INFLUENZA A&B TEST
- INFLUENZA A&B CONTROL KIT

CONFIGURATION

25 Tests

CATALOG NUMBER

190 or 190E (CE marked)
191 or 191E (CE marked)



THE AMERICAS
31 New York Avenue
Framingham, MA 01701-9322
Phone: 800 332 1042
Fax: 800 762 6311
Email: usdiagcustomerservice@genzyme.com

INTERNATIONAL
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF, UK
Phone: +44 (0) 1732 220022
Fax: +44 (0) 1732 220024/5
Email: ukdiagcustomerservice@genzyme.com

